

P2Y4 Polyclonal Antibody

Catalog No: #41302

Package Size: #41302-1 50ul #41302-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	P2Y4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	P2Y4 Polyclonal Antibody detects endogenous levels of P2Y4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human P2Y4.
Target Name	P2Y4
Other Names	P2RY4; NRU; P2Y purinoceptor 4; P2Y4; P2P; Uridine nucleotide receptor; UNR
Accession No.	Swiss-Prot: P51582NCBI Gene ID: 5030
Uniprot	P51582
GeneID	5030;
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using P2Y4 Polyclonal Antibody

Note: This product is for in vitro research use only