

## ARAF antibody

Catalog No: #22653

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	ARAF antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 371 and 524 of ARAF
Target Name	ARAF
Accession No.	Swiss-Prot:P10398Gene ID:369
Uniprot	P10398
GeneID	369;
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details


Predicted MW: 68kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Immunofluorescence: 1:100-1:200

## Images




Sample (30 ug of whole cell lysate)


A: HeLa

7.5% SDS PAGE


Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Cal27 xenograft, using ARAF antibody at 1: 500 dilution.



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed A431, using ARAF antibody at 1: 200 dilution.

## Background

This proto-oncogene belongs to the RAF subfamily of the Ser/Thr protein kinase family, and maybe involved in cell growth and development. [provided by RefSeq]

Note: This product is for in vitro research use only