

## p73 Polyclonal Antibody

Catalog No: #41319

Package Size: #41319-1 50ul #41319-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	p73 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	p73 Polyclonal Antibody detects endogenous levels of p73 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human p73 around the non-phosphorylation site of Y99.
Target Name	p73
Other Names	TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein
Accession No.	Swiss-Prot: O15350NCBI Gene ID: 7161
Uniprot	O15350
GeneID	7161;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of 293T cells using p73 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using p73 Polyclonal Antibody

---

Note: This product is for in vitro research use only