

Pax-1 Polyclonal Antibody

Catalog No: #41322

Package Size: #41322-1 50ul #41322-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pax-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	Pax-1 Polyclonal Antibody detects endogenous levels of Pax-1 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Pax-1.
Target Name	Pax-1
Other Names	PAX1; HUP48; Paired box protein Pax-1; HuP48
Accession No.	Swiss-Prot: P15863NCBI Gene ID: 5075
Uniprot	P15863
GeneID	5075;
SDS-PAGE MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using Pax-1 Polyclonal Antibody

Note: This product is for in vitro research use only