

PBK Polyclonal Antibody

Catalog No: #41325

Package Size: #41325-1 50ul #41325-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PBK Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Mk
Specificity	PBK Polyclonal Antibody detects endogenous levels of PBK protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human PBK around the non-phosphorylation site of T9.
Target Name	PBK
Other Names	PBK; TOPK; Lymphokine-activated killer T-cell-originated protein kinase; Cancer/testis antigen 84; CT84; MAPKK-like protein kinase; Nori-3; PDZ-binding kinase; Spermatogenesis-related protein kinase; SPK; T-LAK cell-originated protein kinas
Accession No.	Swiss-Prot: Q96KB5NCBI Gene ID: 55872
Uniprot	Q96KB5
GeneID	55872;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images

Western Blot analysis of 293T cells using PBK Polyclonal Antibody



Note: This product is for in vitro research use only