

PD2R Polyclonal Antibody

Catalog No: #41326

Package Size: #41326-1 50ul #41326-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PD2R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	PD2R Polyclonal Antibody detects endogenous levels of PD2R protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human PD2R.
Target Name	PD2R
Other Names	PTGDR; Prostaglandin D2 receptor; PGD receptor; PGD2 receptor; Prostanoid DP receptor
Accession No.	Swiss-Prot: Q13258NCBI Gene ID: 5729
Uniprot	Q13258
GeneID	5729;
SDS-PAGE MW	40kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

Images

Western Blot analysis of 293 cells using PD2R Polyclonal Antibody



Note: This product is for in vitro research use only