

PDHA1 Polyclonal Antibody

Catalog No: #41328

Package Size: #41328-1 50ul #41328-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PDHA1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	PDHA1 Polyclonal Antibody detects endogenous levels of PDHA1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human PDHA1.
Target Name	PDHA1
Other Names	PDHA1; PHE1A; Pyruvate dehydrogenase E1 component subunit alpha; somatic form, mitochondrial; PDHE1-A type I
Accession No.	Swiss-Prot: P08559NCBI Gene ID: 5160
Uniprot	P08559
GeneID	5160;
SDS-PAGE MW	43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images

Western Blot analysis of Jurkat cells using PDHA1 Polyclonal Antibody



Note: This product is for in vitro research use only