

PDK1 Polyclonal Antibody

Catalog No: #41329

Package Size: #41329-1 50ul #41329-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PDK1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	PDK1 Polyclonal Antibody detects endogenous levels of PDK1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human PDK1 around the non-phosphorylation site of S241.
Target Name	PDK1
Other Names	PDPK1; PDK1; 3-phosphoinositide-dependent protein kinase 1; hPDK1
Accession No.	Swiss-Prot: O15530NCBI Gene ID: 5170
Uniprot	O15530
GeneID	5170;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.


Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using PDK1 Polyclonal Antibody

Western Blot analysis of HELA COLO 293T cells using PDK1 Polyclonal Antibody



Note: This product is for in vitro research use only