

PDK2 Polyclonal Antibody

Catalog No: #41330

Package Size: #41330-1 50ul #41330-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PDK2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	PDK2 Polyclonal Antibody detects endogenous levels of PDK2 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human PDK2.
Target Name	PDK2
Other Names	PDK2; PDHK2; [Pyruvate dehydrogenase [lipoamide]] kinase isozyme 2; mitochondrial; Pyruvate dehydrogenase kinase isoform 2
Accession No.	Swiss-Prot: Q15119NCBI Gene ID: 5164
Uniprot	Q15119
GeneID	5164;
SDS-PAGE MW	46kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000 IHC 1:100 - 1:300 ELISA: 1:20000 IF 1:50-200

Images



Western Blot analysis of HepG2 cells using PDK2 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

Note: This product is for in vitro research use only