

PEPT1 Polyclonal Antibody

Catalog No: #41331

Package Size: #41331-1 50ul #41331-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PEPT1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	PEPT1 Polyclonal Antibody detects endogenous levels of PEPT1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human PEPT1.
Target Name	PEPT1
Other Names	SLC15A1; PEPT1; Solute carrier family 15 member 1; Intestinal H(+)/peptide cotransporter; Oligopeptide transporter; small intestine isoform; Peptide transporter 1
Accession No.	Swiss-Prot: P46059NCBI Gene ID: 6564
Uniprot	P46059
GeneID	6564;
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of MOUSE-BRAIN MOUSE-HEART MOUSE-LUNG MOUSE-KIDNEY cells using PEPT1 Polyclonal Antibody

Note: This product is for in vitro research use only