

Peroxin 11 β Polyclonal Antibody

Catalog No: #41332

Package Size: #41332-1 50ul #41332-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Peroxin 11 β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	Peroxin 11 β Polyclonal Antibody detects endogenous levels of Peroxin 11 β protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Peroxin 11 β .
Target Name	Peroxin 11 β
Other Names	PEX11B; Peroxisomal membrane protein 11B; Peroxin-11B; Peroxisomal biogenesis factor 11B; Protein PEX11 homolog beta; PEX11-beta
Accession No.	Swiss-Prot: O96011NCBI Gene ID: 8799
Uniprot	O96011
GeneID	8799;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of Jurkat cells using Peroxin 11 β Polyclonal Antibody

Note: This product is for in vitro research use only