

## PGE synthase Polyclonal Antibody

Catalog No: #41334

Package Size: #41334-1 50ul #41334-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PGE synthase Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	PGE synthase Polyclonal Antibody detects endogenous levels of PGE synthase protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human PGE synthase.
Target Name	PGE synthase
Other Names	PTGES; MGST1L1; MPGES1; PGES; PIG12; Prostaglandin E synthase; Microsomal glutathione S-transferase 1-like 1; MGST1-L1; Microsomal prostaglandin E synthase 1; MPGES-1; p53-induced gene 12 protein
Accession No.	Swiss-Prot: O14684NCBI Gene ID: 9536
Uniprot	O14684
GeneID	9536;
SDS-PAGE MW	17kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images

Western Blot analysis of A549 cells using PGE synthase Polyclonal Antibody



Note: This product is for in vitro research use only