

## PHF3 Polyclonal Antibody

Catalog No: #41336

Package Size: #41336-1 50ul #41336-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	PHF3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	PHF3 Polyclonal Antibody detects endogenous levels of PHF3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human PHF3.
Target Name	PHF3
Other Names	PHF3; KIAA0244; PHD finger protein 3
Accession No.	Swiss-Prot: Q92576NCBI Gene ID: 23469
Uniprot	Q92576
GeneID	23469;
SDS-PAGE MW	230kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

## Images



Western Blot analysis of HeLa cells using PHF3 Polyclonal Antibody

---

Note: This product is for in vitro research use only