

PI 3-kinase p85 α / γ Polyclonal Antibody

Catalog No: #41340

Package Size: #41340-1 50ul #41340-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | PI 3-kinase p85 α / γ Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications | WB IHC IF ELISA |
| Species Reactivity | Hu Ms Rt Mk |
| Specificity | PI 3-kinase p85 α / γ Polyclonal Antibody detects endogenous levels of PI 3-kinase p85 α / γ protein. |
| Immunogen Type | peptide |
| Immunogen Description | Synthesized peptide derived from human PI 3-kinase p85/p55 around the non-phosphorylation site of Y458/199. |
| Target Name | PI 3-kinase p85 α / γ |
| Other Names | PIK3R1; GRB1; Phosphatidylinositol 3-kinase regulatory subunit alpha; PI3-kinase regulatory subunit alpha; PI3K regulatory subunit alpha; PtdIns-3-kinase regulatory subunit alpha; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit alph |
| Accession No. | Swiss-Prot: P27986/Q92569NCBI Gene ID: 5295/8503 |
| Uniprot | P27986 |
| GeneID | 5295; |
| SDS-PAGE MW | 85kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat HELA cells using PI 3-kinase p85α/γ Polyclonal Antibody

Note: This product is for in vitro research use only