

PIASx Polyclonal Antibody

Catalog No: #41343

Package Size: #41343-1 50ul #41343-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PIASx Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	PIASx Polyclonal Antibody detects endogenous levels of PIASx protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human PIASx.
Target Name	PIASx
Other Names	PIAS2; PIASX; E3 SUMO-protein ligase PIAS2; Androgen receptor-interacting protein 3; ARIP3; DAB2-interacting protein; DIP; Msx-interacting zinc finger protein; Miz1; PIAS-NY protein; Protein inhibitor of activated STAT x; Protein inhibitor
Accession No.	Swiss-Prot: O75928NCBI Gene ID: 9063
Uniprot	O75928
GeneID	9063;
SDS-PAGE MW	68kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using PIASx Polyclonal Antibody



Western Blot analysis of NIH-3T3 Jurkat K562 cells using PIASx Polyclonal Antibody

Note: This product is for in vitro research use only