

PIG-X Polyclonal Antibody

Catalog No: #41344

Package Size: #41344-1 50ul #41344-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PIG-X Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	PIG-X Polyclonal Antibody detects endogenous levels of PIG-X protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human PIG-X.
Target Name	PIG-X
Other Names	PIGX; Phosphatidylinositol-glycan biosynthesis class X protein; PIG-X
Accession No.	Swiss-Prot: Q8TBF5NCBI Gene ID: 54965
Uniprot	Q8TBF5
GeneID	54965;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using PIG-X Polyclonal Antibody

Note: This product is for in vitro research use only