

PKC ϵ Polyclonal Antibody

Catalog No: #41348

Package Size: #41348-1 50ul #41348-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PKC ϵ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	PKC ϵ Polyclonal Antibody detects endogenous levels of PKC ϵ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human PKC ϵ around the non-phosphorylation site of S729.
Target Name	PKC ϵ
Other Names	PRKCE; PKCE; Protein kinase C epsilon type; nPKC-epsilon
Accession No.	Swiss-Prot: Q02156NCBI Gene ID: 5581
Uniprot	Q02156
GeneID	5581;
SDS-PAGE MW	83kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/5000.

Not yet tested in other applications.

Images

Western Blot analysis of 453 cells using PKC ϵ Polyclonal Antibody

Note: This product is for in vitro research use only