

## ZHX3 Antibody

Catalog No: #47486

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZHX3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms
Specificity	The antibody detects endogenous levels of total ZHX3 ? protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human ZHX3 ?
Target Name	ZHX3
Other Names	TIX1
Accession No.	Swiss-Prot#:Q9H4I2NCBI Gene ID:23051Gene Accssion:NP_055850
Uniprot	Q9H4I2
GeneID	23051;
Concentration	1.3
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

IHC dilution:1: 50-200

## Images



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using 47486(ZHX3 ? Antibody) at dilution 1/100.(Original magnification: 200)

## Background

This gene encodes a member of the zinc fingers and homeoboxes (ZHX) gene family. The encoded protein contains two C<sub>2</sub>H<sub>2</sub>-type zinc fingers and five homeodomains and forms a dimer with itself or with zinc fingers and homeoboxes family member 1. In the nucleus, the dimerized protein interacts with the A subunit of the ubiquitous transcription factor nuclear factor- $\kappa$ B and may function as a transcriptional repressor.

---

Note: This product is for in vitro research use only