

PML Polyclonal Antibody

Catalog No: #41354

Package Size: #41354-1 50ul #41354-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PML Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	PML Polyclonal Antibody detects endogenous levels of PML protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human PML.
Target Name	PML
Other Names	PML; MYL; RNF71; TRIM19; Protein PML; Promyelocytic leukemia protein; RING finger protein 71; Tripartite motif-containing protein 19
Accession No.	Swiss-Prot: P29590NCBI Gene ID: 5371
Uniprot	P29590
GeneID	5371;
SDS-PAGE MW	98kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of A549 cells using PML Polyclonal Antibody



Western Blot analysis of 293T-UV 22RV1 cells using PML Polyclonal Antibody

Note: This product is for in vitro research use only