

POTE-14/22 Polyclonal Antibody

Catalog No: #41356

Package Size: #41356-1 50ul #41356-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	POTE-14/22 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	POTE-14/22 Polyclonal Antibody detects endogenous levels of POTE-14/22 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human POTE-14/22.
Target Name	POTE-14/22
Other Names	POTEH; A26C3; ACTBL1; POTE22; POTE ankyrin domain family member H; ANKRD26-like family C member 3; Prostate; ovary, testis-expressed protein on chromosome 22; POTE-22; POTE G; A26C2; POTE14; POTE ankyrin domain family member G; ANKRD26-like
Accession No.	Swiss-Prot: Q6S545/Q6S5H5NCBI Gene ID: 23784/404785
Uniprot	Q6S545
GeneID	23784;
SDS-PAGE MW	61kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Note: This product is for in vitro research use only