

ABHD17B Antibody

Catalog No: #47537

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABHD17B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms, Rt
Specificity	The antibody detects endogenous levels of total ABHD17B protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human ABHD17B
Target Name	ABHD17B
Other Names	CGI-67; C9orf77; FAM108B1
Accession No.	Swiss-Prot#:Q5VST6NCBI Gene ID:51104Gene Accssion:NP_001020951
Uniprot	Q5VST6
GeneID	51104;
Concentration	0.6
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

IHC dilution:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using 47537(ABHD17B Antibody) at dilution 1/20.(Original magnification: 200)

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only