

## PP32R1 Polyclonal Antibody

Catalog No: #41358

Package Size: #41358-1 50ul #41358-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	PP32R1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	PP32R1 Polyclonal Antibody detects endogenous levels of PP32R1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human PP32R1.
Target Name	PP32R1
Other Names	ANP32C; PP32R1; Acidic leucine-rich nuclear phosphoprotein 32 family member C; Phosphoprotein 32-related protein 1; Tumorigenic protein pp32r1
Accession No.	Swiss-Prot: O43423NCBI Gene ID: 23520
Uniprot	O43423
GeneID	23520;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of HuvEc cells using PP32R1 Polyclonal Antibody

---

Note: This product is for in vitro research use only