

PPAR- $\alpha$  Polyclonal Antibody

Catalog No: #41359

Package Size: #41359-1 50ul #41359-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PPAR- $\alpha$ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA IHC
Species Reactivity	Hu Ms Rt
Specificity	PPAR- $\alpha$ Polyclonal Antibody detects endogenous levels of PPAR- $\alpha$ protein.
Immunogen Description	Synthesized peptide derived from human PPAR- $\alpha$ around the non-phosphorylation site of S21.
Target Name	PPAR- $\alpha$
Other Names	PPARA; NR1C1; PPAR; Peroxisome proliferator-activated receptor alpha; PPAR-alpha; Nuclear receptor subfamily 1 group C member 1
Accession No.	Swiss-Prot: Q07869NCBI Gene ID: 5465
Uniprot	Q07869
GeneID	5465;
SDS-PAGE MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

IHC 1:50-300

## Images

Western Blot analysis of 3T3 cells using PPAR- $\alpha$  Polyclonal Antibody



Note: This product is for in vitro research use only