

PPAR- γ Polyclonal Antibody

Catalog No: #41360

Package Size: #41360-1 50ul #41360-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PPAR- γ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF/ICC
Species Reactivity	Hu Ms Rt
Specificity	PPAR- γ Polyclonal Antibody detects endogenous levels of PPAR- γ protein.
Immunogen Description	Synthesized peptide derived from human PPAR- γ around the non-phosphorylation site of S112.
Target Name	PPAR- γ
Other Names	PPARG; NR1C3; Peroxisome proliferator-activated receptor gamma; PPAR-gamma; Nuclear receptor subfamily 1 group C member 3
Accession No.	Swiss-Prot: P37231NCBI Gene ID: 5468
Uniprot	P37231
GeneID	5468;
SDS-PAGE MW	57kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.


Not yet tested in other applications.

IF: 1:50-200

Images



Western Blot analysis of HuvEc cells using PPAR- γ Polyclonal Antibody



Western Blot analysis of various cells using PPAR- γ Polyclonal Antibody

Note: This product is for in vitro research use only