

PRAF2 Polyclonal Antibody

Catalog No: #41361

Package Size: #41361-1 50ul #41361-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PRAF2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF
Species Reactivity	Hu Ms Rt
Specificity	PRAF2 Polyclonal Antibody detects endogenous levels of PRAF2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human PRAF2.
Target Name	PRAF2
Other Names	PRAF2; JM4; PRA1 family protein 2
Accession No.	Swiss-Prot: O60831NCBI Gene ID: 11230
Uniprot	O60831
GeneID	11230;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of RAT-HEART cells using PRAF2 Polyclonal Antibody



Western Blot analysis of Jurkat NIH-3T3 cells using PRAF2 Polyclonal Antibody

Note: This product is for in vitro research use only