

PRPH Polyclonal Antibody

Catalog No: #41362

Package Size: #41362-1 50ul #41362-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PRPH Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Rt
Specificity	PRPH Polyclonal Antibody detects endogenous levels of PRPH protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human PRPH.
Target Name	PRPH
Other Names	PRPH; NEF4; PRPH1; Peripherin; Neurofilament 4
Accession No.	Swiss-Prot: P41219NCBI Gene ID: 5630
Uniprot	P41219
GeneID	5630;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.


Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using PRPH Polyclonal Antibody

Western Blot analysis of HepG2 cells using PRPH Polyclonal Antibody



Note: This product is for in vitro research use only