

PRX III Polyclonal Antibody

Catalog No: #41363

Package Size: #41363-1 50ul #41363-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | PRX III Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications | WB IHC ELISA |
| Species Reactivity | Hu Ms Rt |
| Specificity | PRX III Polyclonal Antibody detects endogenous levels of PRX III protein. |
| Immunogen Type | peptide |
| Immunogen Description | Synthesized peptide derived from the Internal region of human PRX III. |
| Target Name | PRX III |
| Other Names | PRDX3; AOP1; Thioredoxin-dependent peroxide reductase; mitochondrial; Antioxidant protein 1; AOP-1; HBC189; Peroxiredoxin III; Prx-III; Peroxiredoxin-3; Protein MER5 homolog |
| Accession No. | Swiss-Prot: P30048NCBI Gene ID: 10935 |
| Uniprot | P30048 |
| GeneID | 10935; |
| SDS-PAGE MW | 26kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300. Immunocytochemistry: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images

Western Blot analysis of HeLa cells using PRX III Polyclonal Antibody



Note: This product is for in vitro research use only