

GPR160 Antibody

Catalog No: #47605

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPR160 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GPR160 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human GPR160
Target Name	GPR160
Other Names	GPCR1; GPCR150
Accession No.	Swiss-Prot#:Q9UJ42NCBI Gene ID:26996Gene Accssion:NP_055188
Uniprot	Q9UJ42
GeneID	26996;
Concentration	0.2
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

IHC dilution:1: 30-150

Images



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using 47605(GPR160 Antibody) at dilution 1/20.(Original magnification: 200)

Background

Orphan receptor.

Note: This product is for in vitro research use only