

PTH/PTHrP-R Polyclonal Antibody

Catalog No: #41368

Package Size: #41368-1 50ul #41368-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PTH/PTHrP-R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	PTH/PTHrP-R Polyclonal Antibody detects endogenous levels of PTH/PTHrP-R protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human PTH/PTHrP-R.
Target Name	PTH/PTHrP-R
Other Names	PTH1R; PTHR; PTHR1; Parathyroid hormone/parathyroid hormone-related peptide receptor; PTH/PTHrP type I receptor; PTH/PTHr receptor; Parathyroid hormone 1 receptor; PTH1 receptor
Accession No.	Swiss-Prot: Q03431NCBI Gene ID: 5745
Uniprot	Q03431
GeneID	5745;
SDS-PAGE MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of COLO cells using PTH/PTHrP-R Polyclonal Antibody



Note: This product is for in vitro research use only