

PTP1B Polyclonal Antibody

Catalog No: #41369

Package Size: #41369-1 50ul #41369-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PTP1B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	PTP1B Polyclonal Antibody detects endogenous levels of PTP1B protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human PTP1B around the non-phosphorylation site of S50.
Target Name	PTP1B
Other Names	PTPN1; PTP1B; Tyrosine-protein phosphatase non-receptor type 1; Protein-tyrosine phosphatase 1B; PTP-1B
Accession No.	Swiss-Prot: P18031NCBI Gene ID: 5770
Uniprot	P18031
GeneID	5770;
SDS-PAGE MW	49kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using PTP1B Polyclonal Antibody

Note: This product is for in vitro research use only