

FOXK2 Antibody

Catalog No: #47654



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FOXK2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IHC
Species Reactivity	Hu, Ms
Specificity	The antibody detects endogenous levels of total FOXK2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human FOXK2
Target Name	FOXK2
Other Names	ILF; ILF1; ILF-1
Accession No.	Swiss-Prot#:Q01167NCBI Gene ID:3607Gene Accssion:NP_004505
Uniprot	Q01167
GeneID	3607;
Calculated MW	69 kDa
Concentration	1.4
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN3, 40% Glycerol.
Storage	Store at -20°C


Application Details

WB dilution:1:1000-5000
IHC dilution:1: 50-300

Images



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane: TM4 cell lysate,
Primary antibody:47654(FOXK2 Antibody) at dilution 1/1000,
Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution,
Exposure time: 1 minute



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using 47654(FOXP2 Antibody) at dilution 1/50.(Original magnification: 200)

Background

The protein encoded by this gene contains a fork head DNA binding domain. This protein can bind to the purine-rich motifs of the HIV long terminal repeat (LTR), and to the similar purine-rich motif in the interleukin 2 (IL2) promoter. It may be involved in the regulation of viral and cellular promoter elements.

Note: This product is for in vitro research use only