

PTP ζ Polyclonal Antibody

Catalog No: #41370

Package Size: #41370-1 50ul #41370-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PTP ζ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	PTP ζ Polyclonal Antibody detects endogenous levels of PTP ζ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human PTP ζ .
Target Name	PTP ζ
Other Names	PTPRZ1; HTPZP2; PTPRZ; PTPRZ2; PTPZ; Receptor-type tyrosine-protein phosphatase zeta; R-PTP-zeta; Protein-tyrosine phosphatase receptor type Z polypeptide 1; Protein-tyrosine phosphatase receptor type Z polypeptide 2; R-PTP-zeta-2
Accession No.	Swiss-Prot: P23471NCBI Gene ID: 5803
Uniprot	P23471
GeneID	5803;
SDS-PAGE MW	250kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of A549 cells using PTP ζ Polyclonal Antibody



Western Blot analysis of 293T cells using PTP ζ Polyclonal Antibody

Note: This product is for in vitro research use only