

QM Polyclonal Antibody

Catalog No: #41374



Package Size: #41374-1 50ul #41374-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	QM Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	QM Polyclonal Antibody detects endogenous levels of QM protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human QM.
Target Name	QM
Other Names	RPL10; DXS648E; QM; 60S ribosomal protein L10; Laminin receptor homolog; Protein QM; Tumor suppressor QM
Accession No.	Swiss-Prot: P27635NCBI Gene ID: 6134
Uniprot	P27635
GeneID	6134;
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Note: This product is for in vitro research use only