

Rab 20 Polyclonal Antibody

Catalog No: #41375

Package Size: #41375-1 50ul #41375-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rab 20 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	Rab 20 Polyclonal Antibody detects endogenous levels of Rab 20 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Rab 20.
Target Name	Rab 20
Other Names	RAB20; Ras-related protein Rab-20
Accession No.	Swiss-Prot: Q9NX57NCBI Gene ID: 55647
Uniprot	Q9NX57
GeneID	55647;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of A549 cells using Rab 20 Polyclonal Antibody



Western Blot analysis of Jurkat K562 NIH-3T3 cells using Rab 20 Polyclonal Antibody

Note: This product is for in vitro research use only