

Rab 3 GAP p130 Polyclonal Antibody

Catalog No: #41376

Package Size: #41376-1 50ul #41376-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rab 3 GAP p130 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	Rab 3 GAP p130 Polyclonal Antibody detects endogenous levels of Rab 3 GAP p130 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Rab 3 GAP p130.
Target Name	Rab 3 GAP p130
Other Names	RAB3GAP1; KIAA0066; RAB3GAP; Rab3 GTPase-activating protein catalytic subunit; RAB3 GTPase-activating protein 130 kDa subunit; Rab3-GAP p130; Rab3-GAP
Accession No.	Swiss-Prot: Q15042NCBI Gene ID: 22930
Uniprot	Q15042
GeneID	22930;
SDS-PAGE MW	117kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Rab 3 GAP p130 Polyclonal Antibody

Note: This product is for in vitro research use only