

Rab 34 Polyclonal Antibody

Catalog No: #41377

Package Size: #41377-1 50ul #41377-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rab 34 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	Rab 34 Polyclonal Antibody detects endogenous levels of Rab 34 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Rab 34.
Target Name	Rab 34
Other Names	RAB34; RAB39; RAH; Ras-related protein Rab-34; Ras-related protein Rab-39; Ras-related protein Rah
Accession No.	Swiss-Prot: Q9BZG1NCBI Gene ID: 83871
Uniprot	P0DI83
GeneID	83871;
SDS-PAGE MW	29kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Rab 34 Polyclonal Antibody

Note: This product is for in vitro research use only