

Rab 7L1 Polyclonal Antibody

Catalog No: #41379

Package Size: #41379-1 50ul #41379-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Rab 7L1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Rab 7L1 Polyclonal Antibody detects endogenous levels of Rab 7L1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Rab 7L1.
Target Name	Rab 7L1
Other Names	RAB7L1; Ras-related protein Rab-7L1; Rab-7-like protein 1
Accession No.	Swiss-Prot: O14966NCBI Gene ID: 8934
Uniprot	O14966
GeneID	8934;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HT29 cells using Rab 7L1 Polyclonal Antibody

Note: This product is for in vitro research use only