

## Rad51 Polyclonal Antibody

Catalog No: #41381

Package Size: #41381-1 50ul #41381-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Rad51 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Rad51 Polyclonal Antibody detects endogenous levels of Rad51 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Rad51 around the non-phosphorylation site of Y315.
Target Name	Rad51
Other Names	RAD51; RAD51A; RECA; DNA repair protein RAD51 homolog 1; HsRAD51; hRAD51; RAD51 homolog A
Accession No.	Swiss-Prot: Q06609NCBI Gene ID: 5888
Uniprot	Q06609
GeneID	5888;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.


Not yet tested in other applications.

## Images



Western Blot analysis of HELA-UV cells using Rad51 Polyclonal Antibody

Western Blot analysis of 293T-UV HELA-UV cells using Rad51 Polyclonal Antibody



Note: This product is for in vitro research use only