

## DZIP3 Antibody

Catalog No: #47820

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DZIP3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB, IHC
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of DZIP3 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human DZIP3.
Target Name	DZIP3
Other Names	KIAA0675; E3 ubiquitin-protein ligase DZIP3; DAZ-interacting protein 3; RNA-binding ubiquitin ligase of 138 kDa; hRUL138
Accession No.	Swiss-Prot#:Q86Y13NCBI Gene ID:9666
Uniprot	Q86Y13
GeneID	9666;
Calculated MW	180139KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C


## Application Details

WB (1/500 - 1/1000), IHC (1/100 - 1/200)

## Images



Immunohistochemical analysis of DZIP3 staining in human prostate cancer formalin fixed paraffin embedded tissue section. The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.



Western blot analysis of DZIP3 expression in HEK293T (A), LO2 (B) whole cell lysates.

Note: This product is for in vitro research use only