

## MRGX3 Antibody

Catalog No: #47834

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MRGX3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of MRGX3 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human MRGX3.
Target Name	MRGPRX3
Other Names	MRGX3; SNSR1; SNSR2; Mas-related G-protein coupled receptor member X3; Sensory neuron-specific G-protein coupled receptor 1/2
Accession No.	Swiss-Prot#:Q96LB0NCBI Gene ID:117195
Uniprot	Q96LB0
GeneID	117195;
Calculated MW	36KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/1000)

## Images



Western blot analysis of MRGX3 expression in human brain  
(A) whole cell lysates

Note: This product is for in vitro research use only