

PABPC5 Antibody

Catalog No: #47844

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PABPC5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB, IHC
Species Reactivity	Hu,Rt
Specificity	Recognizes endogenous levels of PABPC5 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human PABPC5.
Target Name	PABPC5
Other Names	PABP5; Polyadenylate-binding protein 5; PABP-5; Poly(A)-binding protein 5
Accession No.	Swiss-Prot#:Q96DU9NCBI Gene ID:140886
Uniprot	Q96DU9
GeneID	140886;
Calculated MW	43,39KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C


Application Details

WB (1/500 - 1/1000), IHC (1/100 - 1/200)

Images



Immunohistochemical analysis of PABPC5 staining in human lung cancer formalin fixed paraffin embedded tissue section. The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.



Western blot analysis of PABPC5 expression in THP1 (A), PC3 (B), HeLa (C), rat brain (D) whole cell lysates.

Note: This product is for in vitro research use only