

CD32c Antibody

Catalog No: #47869

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD32c Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of CD32c protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human CD32c.
Target Name	FCGR2C
Other Names	CD32; FCG2; IGFR2; Low affinity immunoglobulin gamma Fc region receptor II-c; IgG Fc receptor II-c; CDw32; Fc-gamma RII-c; Fc-gamma-RIIc; FcRII-c; CD32
Accession No.	Swiss-Prot#:P31995NCBI Gene ID:9103
Uniprot	P31995
GeneID	9103;
Calculated MW	36KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of CD32c expression in HEK293T (A), mouse lung (B) whole cell lysates.

Note: This product is for in vitro research use only