

FNTB Antibody

Catalog No: #47875

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FNTB Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of FNTB protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human FNTB.
Target Name	FNTB
Other Names	Protein farnesyltransferase subunit beta; FTase-beta; CAAX farnesyltransferase subunit beta; Ras proteins prenyltransferase subunit beta
Accession No.	Swiss-Prot#:P49356NCBI Gene ID:2342
Uniprot	P49356
GeneID	100529261;2342;
Calculated MW	49KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of FNTB expression in zebrafish (A) whole cell lysates

Note: This product is for in vitro research use only