

ADAMTS1 Antibody

Catalog No: #47877

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ADAMTS1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of ADAMTS1 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ADAMTS1.
Target Name	ADAMTS1
Other Names	KIAA1346; METH1; A disintegrin and metalloproteinase with thrombospondin motifs 1; ADAM-TS 1; ADAM-TS1; ADAMTS-1; METH-1
Accession No.	Swiss-Prot#:Q9UHI8NCBI Gene ID:9510
Uniprot	Q9UHI8
GeneID	9510;
Calculated MW	110KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of ADAMTS1 expression in A549 (A), H446 (B) whole cell lysates.

Note: This product is for in vitro research use only