

RFWD3 Polyclonal Antibody

Catalog No: #41395

Package Size: #41395-1 50ul #41395-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	RFWD3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	RFWD3 Polyclonal Antibody detects endogenous levels of RFWD3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human RFWD3.
Target Name	RFWD3
Other Names	RFWD3; RNF201; E3 ubiquitin-protein ligase RFWD3; RING finger and WD repeat domain-containing protein 3; RING finger protein 201
Accession No.	Swiss-Prot: Q6PCD5NCBI Gene ID: 55159
Uniprot	Q6PCD5
GeneID	55159;
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of 823 cells using RFWD3 Polyclonal Antibody

Note: This product is for in vitro research use only