

## GPCR19 Antibody

Catalog No: #47886

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GPCR19 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of GPCR19 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human GPCR19.
Target Name	GPBAR1
Other Names	TGR5; G-protein coupled bile acid receptor 1; G-protein coupled receptor GPCR19; hGPCR19; Membrane-type receptor for bile acids; M-BAR; hBG37; BG37
Accession No.	Swiss-Prot#:Q8TDU6NCBI Gene ID:151306
Uniprot	Q8TDU6
GeneID	151306;
Calculated MW	35KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/1000)

## Images



Western blot analysis of GPCR19 expression in HeLa (A), HCT116 (B), A549 (C) whole cell lysates.

Note: This product is for in vitro research use only