

Ubiquilin 3 Antibody

Catalog No: #47912

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Ubiquilin 3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB, IHC
Species Reactivity	Hu,Rt
Specificity	Recognizes endogenous levels of Ubiquilin 3 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Ubiquilin 3.
Target Name	UBQLN3
Other Names	Ubiquilin-3
Accession No.	Swiss-Prot#:Q9H347NCBI Gene ID:50613
Uniprot	Q9H347
GeneID	50613;
Calculated MW	71KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C


Application Details

WB (1/500 - 1/1000), IHC (1/100 - 1/200)

Images



Immunohistochemical analysis of Ubiquilin 3 staining in human testis formalin fixed paraffin embedded tissue section. The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.



Western blot analysis of Ubiquilin 3 expression in A549 (A) whole cell lysates

Note: This product is for in vitro research use only