

GABRA6 Antibody

Catalog No: #47918

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GABRA6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of GABRA6 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human GABRA6.
Target Name	GABRA6
Other Names	Gamma-aminobutyric acid receptor subunit alpha-6; GABA(A) receptor subunit alpha-6
Accession No.	Swiss-Prot#:Q16445NCBI Gene ID:2559
Uniprot	Q16445
GeneID	2559;
Calculated MW	51KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of GABRA6 expression in Jurkat (A), mouse brain (B) whole cell lysates

Note: This product is for in vitro research use only