

## Ribosomal Protein L17 Polyclonal Antibody

Catalog No: #41401

Package Size: #41401-1 50ul #41401-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ribosomal Protein L17 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	Ribosomal Protein L17 Polyclonal Antibody detects endogenous levels of Ribosomal Protein L17 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ribosomal Protein L17.
Target Name	Ribosomal Protein L17
Other Names	RPL17; 60S ribosomal protein L17; 60S ribosomal protein L23; PD-1
Accession No.	Swiss-Prot: P18621NCBI Gene ID: 6139
Uniprot	P18621
GeneID	6139;
SDS-PAGE MW	24kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of COS7 cells using Ribosomal Protein L17 Polyclonal Antibody

---

Note: This product is for in vitro research use only