

## Ribosomal Protein L18 Polyclonal Antibody

Catalog No: #41402

Package Size: #41402-1 50ul #41402-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ribosomal Protein L18 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Ribosomal Protein L18 Polyclonal Antibody detects endogenous levels of Ribosomal Protein L18 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ribosomal Protein L18.
Target Name	Ribosomal Protein L18
Other Names	RPL18; 60S ribosomal protein L18
Accession No.	Swiss-Prot: Q07020NCBI Gene ID: 6141
Uniprot	Q07020
GeneID	6141;
SDS-PAGE MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

## Images



Western Blot analysis of COLO205 cells using Ribosomal Protein L18 Polyclonal Antibody

---

Note: This product is for in vitro research use only